

## Anti-Socs5 Polyclonal Antibody For Wb, Ihc-P, Elisa - 075159

Item Number: CM0025126

### Introduction

Rabbit polyclonal antibody against human SOCS5, validated for WB, IHC-P, ELISA, with cross-reactivity to human and mouse. Ideal for cytokine signaling and ubiquitination studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SOCS5 Rabbit pAb                                                                                                                      |
| Remarks/Alias    | CIS6; CISH6; Cish5; SOCS-5; SOCS5                                                                                                     |
| Species          | Human                                                                                                                                 |
| GeneID (human)   | 9655                                                                                                                                  |
| GeneID           | 9655                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human SOCS5 (NP_054730.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                                 |
| Application      | WB, IHC-P, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse                                                                                                                          |
| SWISS            | O75159                                                                                                                                |
| Protein Weight   | 61kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |