


Anti-Sod1 Rabbit Monoclonal Antibody For Wb, Elisa - P00441

Item Number: CM0003655



Introduction

Rabbit monoclonal antibody targeting human SOD1 (Superoxide dismutase [Cu-Zn]), validated for WB and ELISA. Ideal for studying oxidative stress, antioxidant mechanisms, and neurodegenerative diseases like ALS.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SOD1 Rabbit mAb
Remarks/Alias	ALS; SOD; ALS1; IPOA; STAHP; hSod1; HEL-S-44; homodimer; SOD1
Species	Human
GeneID (Human)	6647
GeneID	6647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human SOD1 (NP_000445.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00441
Protein Weight	16kDa
Shipping	Ice bag