


Anti-Sod1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P00441

Item Number: CM0007586



Introduction

Anti-SOD1 rabbit monoclonal antibody specifically recognizes human superoxide dismutase [Cu-Zn] (SOD1). Suitable for WB, IHC-P, and ELISA applications in oxidative stress and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	SOD1 Rabbit mAb
Remarks/Alias	ALS; SOD; ALS1; IPOA; STAHP; hSod1; HEL-S-44; homodimer; SOD1
Species	Human
Gene ID (Human)	6647
Gene ID	6647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human SOD1 (NP_000445.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P00441
Protein Weight	16kDa
Shipping	Ice bag