


Anti-Sod1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P00441

Item Number: CM0017307



Introduction

Rabbit polyclonal antibody targeting human SOD1 (Cu-Zn superoxide dismutase). Validated in WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Eliminates toxic radicals and detoxifies hydrogen sulfide. Ideal for oxidative stress and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	SOD1 Rabbit pAb
Remarks/Alias	ALS; SOD; ALS1; IPOA; STAHP; hSod1; HEL-S-44; homodimer; SOD1
Species	Human
GeneID (Human)	6647
GeneID	6647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human SOD1 (NP_000445.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00441
Protein Weight	16kDa
Shipping	Ice bag