


Anti-Sod2 Monoclonal Antibody [Ko Validated] For Wb, If-P, Ihc-P, Elisa - P04179

Item Number: CM0002380



Introduction

High-specificity rabbit monoclonal antibody targeting human SOD2 (MnSOD). KO validated. Applications: WB, IF-P, IHC-P, ELISA. Species reactivity: human, mouse, rat. Ideal for mitochondrial oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SOD2 Rabbit mAb
Remarks/Alias	GC1; IPOB; IPO-B; MNSOD; MVCD6; GClnc1; Mn-SOD; D2
Species	Human
GeneID (human)	6648
GeneID	6648
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOD2 (P04179).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04179
Protein Weight	25kDa
Shipping	Ice bag