


Anti-Sod2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04179

Item Number: CM0020359



Introduction

SOD2 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat superoxide dismutase [Mn], a key mitochondrial antioxidant enzyme. Ideal for oxidative stress and disease research.

[Learn More](#)

Parameter	Value
Product Name	SOD2 Rabbit pAb
Remarks/Alias	GC1; IPOB; IPO-B; MNSOD; MVCD6; GClnc1; Mn-SOD; SOD2
Species	Human
GeneID (Human)	6648
GeneID	6648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-222 of human SOD2 (NP_000627.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P04179
Protein Weight	25kDa
Shipping	Ice bag