

Anti-Soga1 Rabbit Polyclonal Antibody For Wb, Elisa - 094964

Item Number: CM0020361

Introduction

CamelBio's Anti-SOGA1 Rabbit Polyclonal Antibody (WB, ELISA) recognizes human, mouse, and rat SOGA1 (MTCL2), a microtubule cross-linking factor involved in mitosis and autophagy. Ideal for research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SOGA1 Rabbit pAb
Remarks/Alias	SOGA; KIAA0889; C20orf117; SOGA1
Species	Human
GeneID (Human)	140710
GeneID	140710
Immunogen	Recombinant Protein; Recombinant fusion protein containing a sequence corresponding to amino acids 650-750 of human SOGA1 (NP_542194.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94964
Protein Weight	160kDa
Shipping	Ice bag