


Anti-Somatostatin Receptor 2 (Sstr2) Rabbit Monoclonal Antibody For Wb, Elisa - P30874

Item Number: CM0002714



Introduction

Rabbit monoclonal antibody targeting Somatostatin Receptor 2 (SSTR2/SRIF-1). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Suitable for GPCR signaling and neuroendocrine tumor studies.

[Learn More](#)

Parameter	Value
Product Name	Somatostatin Receptor 2 (SSTR2) Rabbit mAb
Remarks/Alias	SST2; Somatostatin Receptor 2 (SSTR2)
Species	Human
Gene ID (Human)	6752
Gene ID	6752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 302-369 of human Somatostatin Receptor 2 (SSTR2) (NP_001041.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30874
Protein Weight	41kDa
Shipping	Ice bag