


Anti-Somatostatin (Sst) Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P61278

Item Number: CM0002387



Introduction

Rabbit monoclonal antibody targeting Somatostatin (SST), a key inhibitor of pituitary hormone secretion. Applications: IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Supports neuroendocrine and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Somatostatin (SST) Rabbit mAb
Remarks/Alias	SMST, SST1, Somatostatin (SST)
Species	Human
GeneID (Human)	6750
GeneID	6750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-116 of human Somatostatin (SST) (NP_001039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61278
Protein Weight	13kDa
Shipping	Ice bag