


Anti-Somatostatin (Sst) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P61278

Item Number: CM0028965



Introduction

Rabbit polyclonal antibody targeting Somatostatin (SST). Suitable for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Observed molecular weight ~13 kDa. Ideal for hormone regulation and neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	Somatostatin (SST) Rabbit pAb
Remarks/Alias	SMST; SST1; Somatostatin (SST)
Species	Human
Gene ID	6750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-116 of human SST (NP_001039.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61278
Protein Weight	13kDa
Shipping	Ice bag