

Anti-Somatotropin (Gh1) Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01241

Item Number: CM0025874



Introduction

Rabbit polyclonal antibody against human growth hormone (somatotropin/GH1). For WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse, rat. Ideal for growth control and IGF1 studies.

[Learn More](#)

Parameter	Value
Product Name	GH1 Rabbit pAb
Remarks/Alias	GH; GHN; GH-N; GHB5; IGHD2; hGH-N; IGHD1A; IGHD1B; GH1
Species	Human
Gene ID (Human)	2688
Gene ID	2688
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human GH1 (NP_000506.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01241
Protein Weight	25kDa
Shipping	Ice bag