


Anti-Son Rabbit Polyclonal Antibody (Wb/Elisa) - P18583

Item Number: CM0023425



Introduction

Rabbit polyclonal antibody against human SON protein (264kDa), validated for Western Blot (WB) and ELISA. SON is an RNA-binding protein involved in mRNA splicing of cell-cycle and DNA-repair transcripts. Ideal for gene expression and splicing research.

[Learn More](#)

Parameter	Value
Product Name	SON Rabbit pAb
Remarks/Alias	SON3; BASS1; DBP-5; NREBP; TOKIMS; C21orf50; SON
Species	Human
GeneID (human)	6651
GeneID	6651
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 850-950 of human SON (NP_001278340.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P18583
Protein Weight	264kDa
Shipping	Ice bag