

Anti-Sonic Hedgehog (Shh) Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - Q15465

Item Number: CM0003907



Introduction

CamelBio's Sonic Hedgehog (Shh) Rabbit mAb is a monoclonal antibody for WB, IF/ICC, IF-P, ELISA. Cross-reacts with human, mouse, rat. Immunogen: synthetic peptide. For developmental biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Sonic Hedgehog (Shh) Rabbit mAb
Remarks/Alias	TPT; HHG1; HLP3; HPE3; SMMCI; ShhNC; TPTPS; MCOPCB5; Sonic Hedgehog (Shh)
Species	Human
Gene ID (Human)	6469
Gene ID	6469
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Sonic Hedgehog (Shh) (Q15465).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15465
Protein Weight	50kDa
Shipping	Ice bag