

Anti-Sonic Hedgehog (Shh) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15465

Item Number: CM0013887



Introduction

Rabbit polyclonal antibody against human Sonic Hedgehog (Shh) protein. Validated for WB, IF/ICC, ELISA. For research on developmental patterning, hedgehog signaling, and morphogen function.

[Learn More](#)

Parameter	Value
Product Name	Sonic Hedgehog (Shh) Rabbit pAb
Remarks/Alias	TPT; HHG1; HLP3; HPE3; SMMCI; ShhNC; TPTPS; MCOPCB5; Sonic Hedgehog (Shh)
Species	Human
GeneID	6469
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 231-330 of human Sonic Hedgehog (Shh) (NP_000184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15465
Protein Weight	50kDa
Shipping	Ice bag