

Anti-Sonic Hedgehog (Shh) Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - Q15465

Item Number: CM0025931



Introduction

Rabbit polyclonal antibody targeting human Sonic Hedgehog (Shh), suitable for WB, IF/ICC, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for developmental signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Sonic Hedgehog (Shh) Rabbit pAb
Remarks/Alias	TPT; HHG1; HLP3; HPE3; SMMCI; ShhNC; TPTPS; MCOPCB5; Sonic Hedgehog (Shh)
Species	Human
Gene ID (Human)	6469
Gene ID	6469
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human Sonic Hedgehog (Shh) (NP_000184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15465
Protein Weight	50kDa
Shipping	Ice bag