

Anti-Sorbitol Dehydrogenase Monoclonal Antibody For Wb, Elisa - Q00796

Item Number: CM0006948



Introduction

Accurate anti-Sorbitol Dehydrogenase monoclonal antibody (SORD) validated for WB and ELISA, cross-reacts with human, mouse, rat. For polyol pathway, diabetic complications, and metabolic studies.

[Learn More](#)

Product Name	Sorbitol Dehydrogenase Rabbit mAb
Remarks/Alias	RDH; SDH; XDH; SORD1; SORDD; HEL-S-95n; Sorbitol Dehydrogenase
Species	Human
GeneID (Human)	6652
GeneID	6652
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human Sorbitol Dehydrogenase (NP_003095.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00796
Protein Weight	38kDa
Shipping	Ice bag