

Anti-Sorbitol Dehydrogenase Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q00796

Item Number: CM0020632



Introduction

Rabbit polyclonal antibody targeting Sorbitol Dehydrogenase (SORD, Q00796). Validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. For polyol pathway and diabetic complication studies.

[Learn More](#)

Parameter	Value
Product Name	Sorbitol Dehydrogenase Rabbit pAb
Remarks/Alias	RDH; SDH; XDH; SORD1; SORDD; HEL-S-95n; Sorbitol Dehydrogenase
Species	Human
Gene ID	6652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 145-300 of human Sorbitol Dehydrogenase (NP_003095.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00796
Protein Weight	38 kDa
Shipping	Ice bag