


Anti-Sorl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92673

Item Number: CM0013835



Introduction

Rabbit polyclonal antibody targeting human SORL1 (Q92673), validated for Western blot and ELISA. Sortilin-related receptor is a multi-ligand sorting receptor involved in APP trafficking and Alzheimer's disease, making this antibody ideal for neuroscience and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SORL1 Rabbit pAb
Remarks/Alias	LR11; LRP9; SORLA; gp250; SorLA-1; C11orf32; SORL1
Species	Human
GeneID (Human)	6653
GeneID	6653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2065-2214 of human SORL1 (NP_003096.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92673
Protein Weight	248kDa
Shipping	Ice bag