


Anti-Sort1 Rabbit Polyclonal Antibody For Wb And Elisa - Q99523

Item Number: CM0027046



Introduction

Rabbit polyclonal antibody targeting human sortilin (SORT1). Validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for studies on protein trafficking, lysosomal sorting, neurotensin signaling, and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	SORT1 Rabbit pAb
Remarks/Alias	NT3; Gp95; NTR3; LDLCQ6; SORT1
Species	Human
GeneID (Human)	6272
GeneID	6272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 785-823 of human SORT1 (NP_002950.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99523
Protein Weight	92kDa
Shipping	Ice bag