


Anti-Sos1 Rabbit Polyclonal Antibody For Wb And Elisa - Q07889

Item Number: CM0026961



Introduction

Rabbit polyclonal antibody targeting human SOS1 (Son of sevenless homolog 1). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Suitable for Ras signaling research.

[Learn More](#)

Parameter	Value
Product Name	SOS1 Rabbit pAb
Remarks/Alias	GF1; HGF; NS4; GGF1; GINGF; SOS-1; SOS1
Species	Human
GeneID (Human)	6654
GeneID	6654
Immunogen	Recombinant protein corresponding to amino acids 1034-1333 of human SOS1 (NP_005624.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07889
Protein Weight	152kDa
Shipping	Ice bag