

Anti-Sost Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bqb4

Item Number: CM0019015



Introduction

Rabbit polyclonal antibody targeting human sclerostin (SOST), validated for Western blot and ELISA. Sclerostin is a negative regulator of bone growth via Wnt signaling inhibition, making this antibody suitable for bone remodeling and osteoporosis research.

[Learn More](#)

Parameter	Value
Product Name	SOST Rabbit pAb
Remarks/Alias	SOST; CDD; DAND6; SOST1; VBCH; sclerostin
Species	Human
Gene ID	50964
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-213 of human SOST (NP_079513.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BQB4
Protein Weight	24kDa
Shipping	Ice bag