

Anti-Sost Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bqb4

Item Number: CM0027470



Introduction

SOST Rabbit pAb is a polyclonal antibody against sclerostin, validated for WB and ELISA. It reacts with human and mouse samples. Ideal for bone growth and Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	SOST Rabbit pAb
Remarks/Alias	CDD; VBCH; DAND6; SOST1; SOST
Species	Human
Gene ID (Human)	50964
Gene ID	50964
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 24-213 of human SOST (NP_079513.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BQB4
Protein Weight	24kDa
Shipping	Ice bag