

## Anti-Sostdc1 Polyclonal Antibody For Wb, Elisa - Q6X4U4

Item Number: CM0022776

### Introduction

Rabbit polyclonal antibody against human SOSTDC1 (USAG-1, ectodin), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. A tool for studying BMP/Wnt signaling and reproductive biology.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SOSTDC1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | DAND7; USAG1; CDA019; ECTODIN; SOSTDC1                                                                                |
| Species          | Human                                                                                                                 |
| Gene ID          | 25928                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 115-206 of human SOSTDC1 (NP_056279.1). |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q6X4U4                                                                                                                |
| Protein Weight   | 23kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |