

Anti-Sox1 Monoclonal Antibody For Wb, Elisa - 000570

Item Number: CM0007492

Introduction

Rabbit monoclonal antibody against human SOX1, validated for WB and ELISA. Cross-reacts with mouse and rat. Manufactured using recombinant fusion protein immunogen (aa 1-100). Ideal for neural differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	SOX1 Rabbit mAb
Remarks/Alias	SOX1; SRY-box 1
Species	Human
GeneID (human)	6656
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SOX1. (NP_005977.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00570
Protein Weight	39kDa
Shipping	Ice bag