

Anti-Sox1 Polyclonal Antibody For Wb, Ihc-P, Elisa - 000570

Item Number: CM0030770

Introduction

Affinity-purified rabbit polyclonal antibody against human SOX1, validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with mouse and rat. Suitable for neural development research.

[Learn More](#)

Parameter	Value
Product Name	SOX1 Rabbit pAb
Remarks/Alias	SOX1; SRY-box 1
Species	Human
Gene ID (Human)	6656
Gene ID	6656
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 292-391 of human SOX1 (NP_005977.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00570
Protein Weight	39kDa
Shipping	Ice bag