


Anti-Sox11 Monoclonal Antibody For Wb And Elisa - P35716

Item Number: CM0006561



Introduction

CamelBio offers this SOX11 Rabbit mAb for WB and ELISA; detects human, mouse, and rat. Suitable for transcription factor research and developmental biology. SWISS: P35716.

[Learn More](#)

Parameter	Value
Product Name	SOX11 Rabbit mAb
Remarks/Alias	CSS9; MRD27; IDDMOH; SOX11
Species	Human
GeneID (human)	6664
GeneID	6664
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX11 (NP_003099.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35716
Protein Weight	47kDa
Shipping	Ice bag