

Anti-Sox13 Polyclonal Antibody For Wb And Elisa - Q9Un79

Item Number: CM0022483

Introduction

Rabbit polyclonal antibody targeting human SOX13, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for transcriptional regulation, diabetes, and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	SOX13 Rabbit pAb
Remarks/Alias	SOX13; ICA12; Sox-13; SRY-box 13
Species	Human
GeneID	9580
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 525-600 of human SOX13 (NP_005677.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UN79
Protein Weight	69kDa
Shipping	Ice bag