

Anti-Sox13 Polyclonal Antibody For Wb, Elisa - Q9Un79

Item Number: CM0017128

Introduction

Rabbit polyclonal antibody targeting human SOX13, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for research on transcription regulation and WNT signaling.

[Learn More](#)

Parameter	Value
Product Name	SOX13 Rabbit pAb
Remarks/Alias	ICA12; Sox-13; SOX13
Species	Human
GeneID (Human)	9580
GeneID	9580
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UN79
Protein Weight	69kDa
Shipping	Ice bag