

Anti-Sox17 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H6I2

Item Number: CM0021765



Introduction

Rabbit polyclonal antibody targeting human SOX17, validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studies on transcription regulation, embryonic development, and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	SOX17 Rabbit pAb
Remarks/Alias	VUR3; SOX17
Species	Human
GeneID (Human)	64321
GeneID	64321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 343-414 of human SOX17 (NP_071899.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6I2
Protein Weight	44 kDa
Shipping	Ice bag