


Anti-Sox2 Monoclonal Antibody For Wb, If-F, Elisa - P48431

Item Number: CM0003553



Introduction

Rabbit monoclonal anti-SOX2 antibody validated for WB, IF-F, and ELISA. Detects human SOX2 and cross-reacts with mouse and rat. Suitable for stem cell pluripotency and neural development studies. P48431.

[Learn More](#)

Parameter	Value
Product Name	SOX2 Rabbit mAb
Remarks/Alias	ANOP3; MCOP53; SOX2
Species	Human
GeneID (Human)	6657
GeneID	6657
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX2 (NP_003097.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48431
Protein Weight	34kDa
Shipping	Ice bag