


Anti-Sox2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P48431

Item Number: CM0017432



Introduction

Rabbit polyclonal antibody against human SOX2 (P48431), validated for WB, IHC-P, IF/ICC, IP, and ELISA in human, mouse, and rat. Ideal for stem cell and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SOX2 Rabbit pAb
Remarks/Alias	ANOP3; MCOPS3; SOX2
Species	Human
Gene ID (Human)	6657
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human SOX2 (NP_003097.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48431
Protein Weight	34kDa
Shipping	Ice bag