


Anti-Sox2 Rabbit Polymab® For Wb, Ip, If, Ihc, Chip, Elisa - P48431

Item Number: CM0002640



Introduction

SOX2 Rabbit PolymAb® (P48431) - validated for WB, IP, IF, IHC-P, ChIP, ELISA across human, mouse, rat. Ideal for stem cell research, neural development, and pluripotency studies.

[Learn More](#)

Parameter	Value
Product Name	SOX2 Rabbit PolymAb®
Remarks/Alias	ANOP3; MCOP53
Species	Human
GeneID (Human)	6657
GeneID	6657
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48431
Protein Weight	34 kDa
Shipping	Ice bag