


Anti-Sox30 Rabbit Polyclonal Antibody For Wb, Elisa - 094993

Item Number: CM0018284



Introduction

High-quality anti-SOX30 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat SOX30 (approx. 82 kDa). Ideal for studies on transcription regulation, spermiogenesis, and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	SOX30 Rabbit pAb
Remarks/Alias	SOX30
Species	Human
Gene ID (Human)	11063
Gene ID	11063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 604-753 of human SOX30 (NP_848511.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94993
Protein Weight	82kDa
Shipping	Ice bag