

Anti-Sox30 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Cgw4

Item Number: CM0015139

Introduction

Rabbit polyclonal antibody targeting SOX30, validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for spermiogenesis, Wnt signaling, transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SOX30 Rabbit pAb
Remarks/Alias	SOX30
Species	Mouse
GeneID (Human)	11063
GeneID	214105
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of mouse SOX30 (NP_775560.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8CGW4
Protein Weight	84kDa
Shipping	Ice bag