

Anti-Sox4 Polyclonal Antibody For Wb, Elisa - Q06945

Item Number: CM0017446

Introduction

SOX4 rabbit polyclonal antibody (pAb) for Western blot and ELISA. Cross-reacts with mouse. Targets human transcription factor SOX-4 (47kDa), involved in T-cell development and myogenesis. Ideal for immunological and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SOX4 Rabbit pAb
Remarks/Alias	CSS10; EVI16; SOX4
Species	Human
GeneID (Human)	6659
GeneID	6659
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX4 (NP_003098.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q06945
Protein Weight	47kDa
Shipping	Ice bag