

Anti-Sox5 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35711

Item Number: CM0026773

Introduction

Rabbit polyclonal antibody targeting human SOX5, suitable for western blot, immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for chondrogenesis and transcription research.

[Learn More](#)

Parameter	Value
Product Name	SOX5 Rabbit pAb
Remarks/Alias	L-SOX5; LAMSHF; L-SOX5B; L-SOX5F; SOX5
Species	Human
GeneID	6660
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-120 of human SOX5 (NP_694534.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35711
Protein Weight	84kDa
Shipping	Ice bag