

Anti-Sox5 Rabbit Monoclonal Antibody For Wb, Elisa - P35711

Item Number: CM0006446

Introduction

This SOX5 Rabbit Monoclonal Antibody is suitable for Western blot and ELISA applications. It cross-reacts with mouse and rat, enabling studies of SOX5-mediated chondrocyte differentiation and cartilage development.

[Learn More](#)

Parameter	Value
Product Name	SOX5 Rabbit mAb
Remarks/Alias	L-SOX5; LAMSHF; L-SOX5B; L-SOX5F; SOX5
Species	Human
Gene ID	6660
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX5. (NP_008871.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35711
Protein Weight	84 kDa
Shipping	Ice bag