

## Anti-Sox6 Polyclonal Antibody For If/Icc And Elisa - P35712

Item Number: CM0023685

### Introduction

High-quality rabbit polyclonal antibody against human SOX6 transcription factor, validated for immunofluorescence (IF/ICC) and ELISA. Ideal for studying chondrogenesis, neurogenesis, and developmental processes.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | SOX6 Rabbit pAb                                                                                                 |
| Remarks/Alias    | SOXD; HSSOX6; TOLCAS; SOX6                                                                                      |
| Species          | Human                                                                                                           |
| Gene ID          | 55553                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human PML(NP_001139291.2) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human                                                                                                           |
| SWISS            | P35712                                                                                                          |
| Protein Weight   | 92kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |