

Anti-Sox6 Rabbit Polyclonal Antibody For Wb, Elisa - P35712

Item Number: CM0029111

Introduction

SOX6 rabbit polyclonal antibody (pAb) validated for WB and ELISA applications. Detects human, mouse, and rat SOX6 protein (92kDa), a key transcription factor in development. Use for neurogenesis and chondrogenesis research.

[Learn More](#)

Parameter	Value
Product Name	SOX6 Rabbit pAb
Remarks/Alias	SOXD; HSSOX6; TOLCAS; SOX6
Species	Human
GeneID (Human)	55553
GeneID	55553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-74 of human SOX6 (NP_001139291.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35712
Protein Weight	92kDa
Shipping	Ice bag