

## Anti-Sox9 Rabbit Monoclonal Antibody For Wb - P48436

Item Number: CM0003604



### Introduction

High-quality rabbit monoclonal antibody against human SOX9 transcription factor, validated for WB, IP, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat; ideal for chondrogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                         |
|------------------|---------------------------------------------------------------------------------------------------------------|
| Product Name     | SOX9 Rabbit mAb                                                                                               |
| Remarks/Alias    | CMD1; SRA1; CMPD1; SRXX2; SRXY10; SOX9                                                                        |
| Species          | Human                                                                                                         |
| Gene ID (Human)  | 6662                                                                                                          |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 169-300 of human SOX9 (P48436). |
| Source           | Rabbit                                                                                                        |
| Category         | Monoclonal Antibodies                                                                                         |
| Application      | WB, IP, IF/ICC, IF-P, IHC-P, ELISA                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                             |
| SWISS            | P48436                                                                                                        |
| Protein Weight   | 56kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                       |