

Anti-Sp1 Rabbit Polyclonal Antibody For Ihc-P, Elisa - P08047

Item Number: CM0030198



Introduction

SP1 Rabbit Polyclonal Antibody validated for IHC-P and ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for studying transcription regulation, chromatin remodeling, and cancer biology. High-quality tool for research.

[Learn More](#)

Parameter	Value
Product Name	SP1 Rabbit pAb
Remarks/Alias	SP1
Species	Human
Gene ID (human)	6667
Gene ID	6667
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 581-680 of human SP1 (NP_612482.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08047
Protein Weight	81kDa
Shipping	Ice bag