

Anti-Sp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08047

Item Number: CM0013912



Introduction

Rabbit polyclonal antibody targeting human SP1 transcription factor. Validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat SP1. Suitable for research on cell growth, apoptosis, and chromatin remodeling.

[Learn More](#)

Parameter	Value
Product Name	SP1 Rabbit pAb
Remarks/Alias	SP1
Species	Human
Gene ID (Human)	6667
Gene ID	6667
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 420-620 of human SP1 (NP_612482.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08047
Protein Weight	81kDa
Shipping	Ice bag