


Anti-Sp100 Polyclonal Antibody For Wb And Elisa - P23497

Item Number: CM0026366



Introduction

Rabbit polyclonal SP100 antibody for Western blot and ELISA, reactive with human and mouse. Suitable for investigating PML nuclear bodies, transcriptional coactivation with ETS1, and apoptosis regulation.

[Learn More](#)

Parameter	Value
Product Name	SP100 Rabbit pAb
Remarks/Alias	lysp100b; SP100
Species	Human
GeneID (human)	6672
GeneID	6672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 181-480 of human SP100 (NP_001193631.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23497
Protein Weight	100kDa
Shipping	Ice bag