

Anti-Sp110 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hb58

Item Number: CM0030156

Introduction

Rabbit polyclonal antibody targeting human SP110 (IPR1, IFI41, IFI75). Validated for Western blotting (WB) and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for transcription regulation and autoimmune disease research.

[Learn More](#)

Parameter	Value
Product Name	SP110 Rabbit pAb
Remarks/Alias	IPR1; VODI; IFI41; IFI75; SP110
Species	Human
Gene ID (Human)	3431
Immunogen	Recombinant protein corresponding to amino acids 1-290 of human SP110 (NP_004501.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB58
Protein Weight	78kDa
Shipping	Ice bag