

Anti-Sp2 Polyclonal Antibody For Wb, Elisa - Q02086

Item Number: CM0019766

Introduction

Rabbit polyclonal antibody targeting human SP2 transcription factor (~65 kDa). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	SP2 Rabbit pAb
Species	Human
Gene ID (Human)	6668
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-331 of human SP2 (NP_003101.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q02086
Protein Weight	65kDa
Shipping	Ice bag