

Anti-Sp3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q02447

Item Number: CM0015644

Introduction

Rabbit polyclonal antibody recognizing human, mouse, rat SP3 (Transcription factor Sp3). Validated for WB, IHC-P, ELISA. Suitable for studying transcriptional regulation, cell-cycle control, and gene expression dynamics.

[Learn More](#)

Parameter	Value
Product Name	SP3 Rabbit pAb
Remarks/Alias	SPR2; SP3
Species	Human
GenelD (Human)	6670
GenelD	6670
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02447
Protein Weight	82kDa
Shipping	Ice bag