

Anti-Sp7 (Osterix) Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tdd2

Item Number: CM0019587

Introduction

Anti-SP7 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Recognizes human, mouse, and rat SP7/osterix, a key transcription factor for osteoblast differentiation. Ideal for bone biology research.

[Learn More](#)

Parameter	Value
Product Name	SP7 Rabbit pAb
Remarks/Alias	OSX; OI11; OI12; osterix; SP7
Species	Human
GeneID (Human)	121340
GeneID	121340
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 333-376 of human SP7 (NP_001166938.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDD2
Protein Weight	45kDa
Shipping	Ice bag