

Anti-Spag11A Polyclonal Antibody For Wb, Elisa - Q6Pda7

Item Number: CM0020241

Introduction

Rabbit polyclonal antibody against human SPAG11A, validated for WB and ELISA with cross-reactivity to mouse and rat. Targets a secreted antimicrobial defensin essential for sperm maturation and male reproductive tract defense.

[Learn More](#)

Parameter	Value
Product Name	SPAG11A Rabbit pAb
Remarks/Alias	HE2; EDDM2A; SPAG11A
Species	Human
GeneID (human)	653423
GeneID	653423
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-50 of human SPAG11A (NP_478108.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6PDA7
Protein Weight	14kDa
Shipping	Ice bag