

Anti-Spag5 Polyclonal Antibody For If/Icc, Elisa - Q96R06

Item Number: CM0027647

Introduction

High-quality rabbit polyclonal antibody against human SPAG5 (Astrin, MAP126) for IF/ICC and ELISA. Recombinant immunogen (aa 1-300) ensures specificity for studying mitotic spindle dynamics and chromosome segregation. - Q96R06

[Learn More](#)

Parameter	Value
Product Name	SPAG5 Rabbit pAb
Remarks/Alias	MAP126; DEEPEST; hMAP126; SPAG5
Species	Human
GeneID (Human)	10615
GeneID	10615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SPAG5 (NP_006452.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96R06
Protein Weight	134kDa
Shipping	Ice bag