

Anti-Spag6 Rabbit Polyclonal Antibody For Wb And Elisa - O75602

Item Number: CM0027113

Introduction

Rabbit polyclonal antibody targeting human SPAG6, a sperm-associated antigen essential for flagellar motility. Validated for Western blot and ELISA applications. Suitable for reproductive biology and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	SPAG6 Rabbit pAb
Remarks/Alias	pf16; CT141; FAP194; CFAP194; Repro-SA-1; SPAG6
Species	Human
GeneID (Human)	9576
GeneID	9576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human SPAG6 (NP_758442.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75602
Protein Weight	55kDa
Shipping	Ice bag